

Infectious Disease Testing

State-of-the-Art Diagnosis

CoreLab DX

ELISA


HIV
HIV Ag/Ab
ELISA

Toxo
IgM μ -capture
ELISA

Hepatitis C
anti-HCV
ELISA

Rubella
IgM μ -capture
ELISA

Syphilis
Screen
ELISA



Chlamydia
TRACHOMATIS
IgG, IgA
ELISA

Human

Diagnostics Worldwide

Infectious Disease ELISAs

Superior Quality for Manual or Automated Testing

Perfect for blood bank screening

Product		REF	Sensitivity	Specificity	Sample
HIV Ag/Ab	96 tests	51215	100 %	99,6 %	70 µl serum or plasma
anti-HIV-1/2	96 tests	51205C	100 %	99,6 %	30 µl serum or plasma
anti-HCV	96 tests	51275	100 %	99,7 %	100 µl serum or plasma
HBsAg	96 tests	51048	100 %	99,6 %	50 µl serum or plasma
SYPHILIS SCREEN	96 tests	51005	100 %	100 %	50 µl serum or plasma

Reliable infection serology for prenatal care and routine testing

ToRCH		REF	Other Parameters		REF
CMV (Cytomegalo Virus)	IgG	51203	EBV (Epstein-Barr Virus)	IgG	51204
	IgM	51103		IgM	51104
HSV 1 (Herpes Simplex Virus 1)	IgG	51216	H. PYLORI	IgG	51222
HSV 2 (Herpes Simplex Virus 2)	IgG	51226		IgA	51322
HSV (Herpes Simplex Virus)	IgM	51126	MEASLES VIRUS	IgG	51206
	RUBELLA VIRUS	IgG		IgM	51206
Toxo (Toxoplasma Gondii)	IgG	51209	MUMPS VIRUS	IgG	51207
	IgM µ-capture	51119		IgM	51107
Hepatitis *	Antibody test	51255	VZV (Varicella-Zoster Virus)	IgG	51210
				IgA	51110
anti-HBs	Antibody test	51260	CHLAMYDIA TRACHOMATIS	IgG	51228
				IgA	51328
Tropical Diseases	DENGUE	IgG		IgG	51240
				IgM	51140

Automation for all laboratories

- > The choice of three different instruments allows flexible automation for your laboratory.



* Non CE marked products

981146/2017-10 © 2017



HUMAN Gesellschaft für Biochemica und Diagnostica mbH
 Max-Planck-Ring 21 · 65205 Wiesbaden · Germany
 Tel. +49 6122-9988-0 · Fax +49 6122-9988-100 · e-mail: human@human.de · www.human.de

Human

Diagnostics Worldwide