



TROLLEYS

DRESSING TROLLEYS

Two detachable pressed trays with dimensions 600 x 400 x 60 mm. With rounded corners and polished non-sharp edges.
Two drawers with handle.
Useful area of the drawer: 220 x 340 x 80 mm
Handle to facilitate movements.
Four ø80 mm castors with bumpers, two with brakes.
Made of stainless steel 18/10.

detachable trays
600 x 400 x 60 mm

H•16



OPTIONAL ACCESSORIES

Stainless steel waste bucket with lid. Capacity: 7 litres.
ABS plastic swing bin. Capacity: 9 litres.
Detachable stainless steel wire basket with dimensions 420 x 220 x 150 mm

stainless steel waste
bucket

wire basket



H•16 ∴ 710 x 400 x 830 mm (without accessories) ∴ 15 kg

DRESSING TROLLEYS

Two detachable pressed trays with dimensions 600 x 400 x 60 mm. (capacity: 6 litres). With rounded corners and polished non-sharp edges.

Handle to facilitate movements.

Four ø80 mm castors with bumpers, two with brakes.

Made of stainless steel 18/10.

detachable trays
600 x 400 x 60 mm

H·17



OPTIONAL ACCESSORIES


Stainless steel waste bucket with lid. Capacity: 7 litres.

ABS plastic swing bin. Capacity: 9 litres.

Detachable stainless steel wire basket with dimensions 420 x 220 x 150 mm

wire basket

ABS swing bin



H·17 .. 710 x 400 x 830 mm (without accessories) .. 11 kg



TROLLEYS

DRESSING TROLLEYS

Upper shelf with protection flange (25 mm) on three sides.
With rounded corners and polished non-sharp edges.
Lower flat shelf.
Handle to facilitate movements.
Four $\varnothing 80$ mm castors with bumpers, two with brakes.
Made of stainless steel 18/10.

H-214



model with one drawer
useful area of the drawer: 400 x 340 x 80 mm

H-215



H-216

model with two drawers
useful area of the drawer: 220 x 340 x 80 mm



H-214 .. 710 x 400 x 830 mm .. 8 kg

H-215 .. 710 x 400 x 830 mm .. 11 kg

H-216 .. 710 x 400 x 830 mm .. 12 kg

OPTIONAL ACCESSORIES

Available for models H-214, H-215 and H-216.

Slide removable basket of 420 x 220 x 150 mm, made of stainless steel 18/10.

Stainless steel waste bucket with lid. Capacity: 7 litres.

ABS plastic swing bin. Capacity: 9 litres.

slide removable basket



stainless steel waste bucket



ABS swing bin



DRESSING TROLLEY

Shelves with a 20 mm high perimetrical protective flange.

Working surface of each shelf: 740 x 480 mm

Space between shelves: 750 mm

Rounded and polished corners without sharp edges.

Two front drawers with handle.

Working space of the drawer: 310 x 520 x 120 mm

Double handle for easy movement.

Ø100 mm swivel castors with bumpers; two with brake.

Manufactured in stainless steel, quality 18/10 (AISI 304).

H-291



H-291 .. 830 x 580 x 940 mm .. 22 kg



TROLLEYS

DRESSING TROLLEYS

Shelves with tray shape with dimensions 680 x 420 mm
Rounded corners with no sharp edges.
Perimetrical rim 20 mm height to avoid falls and liquid spills.
Handle to facilitate movements.
Four ø80 mm castors with bumpers, two with brakes.
Made of stainless steel 18/10.

rounded corners

H-225-2

model with two shelves
space between shelves: 500 mm



H-225-3

model with three shelves
space between shelves: 250 mm



H-225-2 .. 750 x 490 x 830 mm .. 8 kg

H-225-3 .. 750 x 490 x 830 mm .. 11 kg

ergonomic handle

model with one drawer
useful area of the drawer: 600 x 400 x 100 mm

H·226-1



H·226-2

model with two drawers
useful area of the drawer: 270 x 400 x 100 mm



H·226-1 .. 750 x 490 x 830 mm ... 12 kg

H·226-2 .. 730 x 420 x 830 mm ... 13 kg