



# AUXILIARY FURNITURE

## WARD SCREENS

Removable and washable plastic cloth.

White is the standard colour but please ask us for other available colours.

Folding panels.

ø50 mm swivel castors.

Structure made of enamelled, chromed or stainless steel.

two panels

H·62



three panels

H·63



H·62 ∴ 1400 - 700 x 400 x 1800 mm ∴ 7 kg

H·63 ∴ 2100 - 700 x 400 x 1800 mm ∴ 10 kg

## WARD SCREENS

Removable and washable plastic cloth.

White is the standard colour but please ask us for other available colours available.

Folding panels.

ø50 mm swivel castors.

Structure made of enamelled, chromed or stainless steel.

four panels



**H·64** ∴ 2800 - 700 x 400 x 1800 mm ∴ 13 kg



# AUXILIARY FURNITURE

## RIGID WARD SCREENS

Rigid panels made of phenolic resin, characterized by its high resistance to impacts, liquids and disinfectants; indestructible and easy to clean.

White is the standard colour but please ask us for other available colours.

Folding panels.

ø65 mm swivel castors.

Structure made of enamelled steel (colour to choose)

two panels

H·262



three panels

H·263



H·262 ∴ 1400 - 700 x 400 x 1800 mm ∴ 17 kg

H·263 ∴ 2100 - 700 x 400 x 1800 mm ∴ 22 kg

## RIGID WARD SCREENS

Rigid panels made of phenolic resin, characterized by its high resistance to impacts, liquids and disinfectants; indestructible and easy to clean.

White is the standard colour but please ask us for other available colours.

Folding panels.

Ø65 mm swivel castors.

Structure made of enamelled steel (colour to choose)

three panels



long-lasting



light weight



resistant



customizable



antistatic



antibacterial



humidity  
resistant



latex-free