



AUXILIARY FURNITURE

I.V. STANDS

I.V. poles with 2 or 4 hooks.

Height adjustable from 1000 to 1950 mm, by manual knob.

ø50 mm double swivel castors.

Consult other castors available.

Structure made of stainless steel 18/10.

Five-leg base with reinforcing ribs.

Refer to next table for base dimensions and materials.

code	diameter	base material	code	diameter	base material
H-161	ø600	stainless steel	H-165	ø500	aluminium
H-162	ø500	stainless steel	H-166	ø600	nylon
H-163	ø500	stainless steel	H-167	ø600	aluminium
H-164	ø500	nylon	H-169	ø600	enamelled steel



ANTI-OVERTURN IV STAND

IV pole with 4 hangers, with load capacity per hanger: 5 kg
 Height adjustable from 1000 to 1950 mm, by manual knob.
 ø25 mm rod made of stainless steel
 ø600 mm anti-overturn base, made of enameled steel, with five ø75 mm swivel castors, two of them with brake with protective bumpers.
 Structure made of stainless steel, quality 18/10.

H-169

swivel support with 4 hangers

H-161

HYDRAULIC IV STAND

Serum holder specially designed for washing and surgical irrigation.
 ø60 mm pole with with swivel support with 4 metal hangers
 Load capacity per hanger: 7 kg
 Height adjustable from 1500 to 2350 mm, by hydraulic pump activated by pedal
 ø600 mm anti-overturn base, made of enameled steel, with five ø50 mm swivel castors, two of them with brake, with protective bumpers.
 Structure made of stainless steel, quality 18/10.

other castors available

Antistatic-conductive
 ø50 mm castor



Standard ø80 mm castor



Antistatic-conductive
 ø80 mm castor

