



WHEELCHAIRS

WHEELCHAIRS

Wheelchair characterized by its great strength and functionality specifically designed for transporting patients in emergency and hospitalization departments.

The structure made of steel pipe is fixed and joined by robust welding, with epoxy enamel finish in standard orange or custom color.

It includes, as standard, an oxygen bottle holder, an IV pole rod holder and an ergonomic and retractable footrest.

The seat and back are upholstered in heavy-duty vinyl fabric over a high-density foam. Upholstery colour is up to choice.

Polyurethane armrests, can be folded down, but cannot be removed to avoid loss.

It includes 2 handles for easy movement, two ø600 mm rear wheels and two ø150 mm front wheels with brake.

All components are easily replaceable.

Safe working load: 150 Kg

· Working width: 500 mm

· Backrest: 420 x 480 mm

· Seat: 470 x 480 mm

· Footrest: 240 x 370 mm

I.V. pole



H-515



ø600 mm rear wheels

oxygen bottle holder



H-515 .. 670 x 880 x 990 mm .. 24 kg

WHEELCHAIRS

Same specifications as model H-515

It includes 2 handles for easy movement, two ø400 mm rear wheels and two ø150 mm front wheels with brake.

Safe working load: 150 Kg

- Working width: 500 mm
- Backrest: 420 x 480 mm
- Seat: 470 x 480 mm
- Footrest: 240 x 370 mm

H-514

ø400 mm rear wheels

ergonomic footrest



WHEELCHAIR SPECIAL MODEL FOR OBESE PATIENTS

Wheelchair characterized by its great strength and functionality specifically designed for transporting obese patients.

The structure made of steel pipe is fixed and joined by robust welding, with epoxy enamel finish in standard violet or custom color.

Safe working load: 300 Kg

Same specifications as model H-515

- Working width: 620 mm
- Backrest: 420 x 580 mm
- Seat: 470 x 580 mm
- Footrest: 240 x 370 mm

wheelchair for obese patients

H-517



H-514 .. 670 x 880 x 990 mm .. 24 kg

H-517 .. 830 x 880 x 990 mm .. 36 kg