



CABINETS

CABINETS WITH CUPBOARD

Upper part: cabinet with double glass door with keylock and two adjustable glass shelves.

Lower part: metallic cupboard with keylock and one fixed metallic shelf.

Made of enamelled or stainless steel 18/10.

H·77



stainless steel



enamelled steel

H·77 ∴ 600 x 300 x 1400 mm ∴ 33 kg

CABINETS WITH CUPBOARD

Upper part: cabinet with double glass door with keylock and two adjustable glass shelves.

Lower part: metallic cupboard with keylock and one fixed metallic shelf.

Made of enamelled or stainless steel 18/10.



stainless steel



enamelled steel

H-78

H-78 .. 1000 x 400 x 1600 mm .. 60 kg



CABINETS

CABINET WITH DRAWERS

Upper part: cabinet with sliding glass doors with keylock and three adjustable glass shelves inside.

Lower part: three drawers.

Made of stainless steel 18/10



H-79 ∴ 650 x 450 x 1600 mm ∴ 65 kg

CABINETS

Glass door(s) with keylock.
 Four adjustable glass shelves.
 Space between shelves: 250 mm
 Made of enamelled steel.
 Base in stainless steel.

H·81

single door



H·82

double door



H·81 .. 600 x 350 x 1500 mm .. 35 kg

H·82 .. 900 x 350 x 1500 mm .. 48 kg



CABINETS

CABINETS

Sliding glass doors with keylock.
Four adjustable glass shelves.
Made of enamelled or stainless steel.
Base in stainless steel.

H-83

enamelled steel



stainless steel

H-83 .. 900 x 350 x 1600 mm .. 65 kg

CABINETS WITH CUPBOARD

Upper part: cabinet with double glass door with keylock and three adjustable glass shelves inside.

Lower part: metallic cupboard with keylock and a fixed metallic shelf.

Made of stainless steel 18/10.

H·85

sliding glass doors



folding doors

model with two drawers and lower cupboard

H·84



H·84 ∴ 900 x 400 x 1800 mm ∴ 90 kg

H·85 ∴ 900 x 400 x 1600 mm ∴ 80 kg



NEWBORN CRIBS

Detachable anatomical basket in transparent thermo-moulded methacrylate.

ø100 mm castors with bumpers, two of them with brake.

Handle for easy movement.

Structure made of stainless steel 18/10

Polyurethane mattress with sanitary cover is included

.. Mattress dimensions: 670 x 360 x 50 mm

H·11

maximum support of the basket
to ensure baby's safety



H·11TR



Trendelenburg and Reverse Trendelenburg
positions (15°) actuated by gas spring

H·11 .. 900 x 480 x 900 mm .. 14 kg

H·11TR .. 900 x 480 x 900 mm .. 16 kg

CHILDREN CRIB

Reinforced frame.

Laterals and front sides made of transparent methacrylate.

Bed with articulated back (60°) and leg (20°) board by crank.

Protected panoramic panels.

Adjustable folding rails.

I.V. pole included to be placed at any side.

Bed made of compact resin.

Structure made of enamelled or stainless steel.

Ø 125 mm castors with bumpers, all of them with brakes

Polyurethane mattress with sanitary cover is included

· Bed dimensions: 1100 x 550 mm

· Mattress dimensions: 1090 x 540 x 80 mm

protected panoramic panels

maximum safety side rails

H·13

Ø125 mm castors with brake



H·13 .. 1250 x 750 x 1400 mm .. 50 kg