

HumaLyte Plus

Ion-Selective Electrolyte Analyzer

- > Na⁺, K⁺, Cl⁻, iCa⁺⁺, pH*
- > Direct ISE technology
- > Touchscreen
- > Autosampler

CLINICAL CHEMISTRY

CoreLab DX



Human

Diagnostics Worldwide

HumaLyte Plus

Flexibility at it's Best

- > Suitable for all kinds of laboratories
- > Easy handling and cost-efficient testing of ISE parameters



Cl⁻ K⁺ Na⁺

REF

HumaLyte Plus³ 17470/10

ISE measurement for

- > Na⁺, K⁺, Cl⁻

One Concept – Two Instruments

- > Analysis time < 60 sec per sample
- > For testing serum, whole blood (heparin), plasma (heparin) or diluted urine
- > Sample volume 150 µl
- > Refillable electrodes
- > Automatic or manual calibration
- > Optional autosampler for 20 samples and STAT
- > LCD touchscreen
- > Quality control statistics and database
- > Internal thermal printer
- > Memory for 200 samples, QC
- > Automatic washing and rinsing



pH* iCa⁺⁺ Cl⁻ K⁺ Na⁺

REF

HumaLyte Plus⁵ 17470/20

ISE measurement for

- > Na⁺, K⁺, Cl⁻, iCa⁺⁺, pH*

System Reagents

REF

Reagent Pack HumaLyte Plus³ 17470/82

Reagent Pack HumaLyte Plus⁵ 17470/83

- > The snap-in replacement technology of its self-contained reagent pack with enclosed waste cartridge ensures reduced exposure to biohazards



Accessories

REF

HumaLyte Autosampler 17470/30

- > Walk-away-capability on demand
- > Easy add-on
- > 20 sample cup positions + 2 STAT

* Used for correction of iCa results only!



Movie

Human

Diagnostics Worldwide



HUMAN Gesellschaft für Biochemica und Diagnostica mbH
Max-Planck-Ring 21 · 65205 Wiesbaden · Germany
Tel. +49 6122-9988-0 · Fax +49 6122-9988-100 · e-mail: human@human.de · www.human.de